

A High-level, Rapid action and Effective disinfectant for:



Surface Disinfection and Aerial Fumigation



Equipment Disinfection

Composition

Each 100gms contains,
Ethylenedioxy) dimethanol.....14.0 gms
Glutaraldehyde.....5.0 gms
Didecyl Dimethyl ammonium chloride.....8.0 gms
Anti-septic Cleansors
Corrosion Inhibitors

- BactiSteril^{PLUS} is a broad-spectrum disinfectant for critical hygiene needs.
- Effective against bacteria, viruses, fungi, and spores.
- Ideal for use in hospitals, OTs, ICUs, laboratories, and cleanrooms.
- Provides thorough decontamination for surfaces and equipment disinfection.
- Aldehyde-based formula with rapid and long-lasting protection.
- Low odor profile for user comfort in closed or populated environments.
- Fast-acting formulation begins microbial inactivation within minutes of application.



Broad-Spectrum

Rapid Action

Effective & Reliable

Long Lasting

DIRECTIONS FOR USE

Optimal Disinfection: Recommended dosage
Add 25ml of BactiSteril^{PLUS} to 5L of water.

OT Disinfection: Add 50ml solution of
BactiSteril^{PLUS} to 5L of water and apply on the
surface or use in the fogger to fumigate.

FIRST AID MEASURES

EYE CONTACT <ul style="list-style-type: none">• Flush eyes with cold water for 15 min.• Seek medical attention immediately.	SERIOUS SKIN CONTACT <ul style="list-style-type: none">• Flush skin with cold water immediately.• Clean with disinfectant soap and apply. antibacterial cream• Seek medical attention.
SERIOUS INHALATION: <ul style="list-style-type: none">• If inhaled - move to fresh air.• If not breathing - give artificial respiration.• Seek medical attention immediately.	INGESTION: <ul style="list-style-type: none">• Do not induce vomiting• Do not give anything by mouth• Seek medical attention immediately.


HANDLING MEASURES & STORAGE

Precautions: Keep away from heat, sources of ignition. Empty containers pose a fire risk, evaporate residue under fume hood. Do not ingest, breathe.

Wear suitable protective clothing. In case of insufficient ventilation, wear suitable respiratory equipments. If ingested seek medical advice. Keep away from oxidizing agents.

Avoid contact with skin and eyes; in case of contact, rinse immediately with plenty of water.

APPLICATIONS



Healthcare



Pharma Labs



Food & Beverage Industries



Research Labs



APPIDI HEALTHCARE, A Division of Appidi Technologies Pvt. Ltd.

Survey no 123, Jeedimetla village,

Quthbullapur Road,

Hyderabad-500055, India.

+91 91546 21871



+91 9154579984



healthcare@appiditech.com



www.appidihealthcare.com